
ForaCare™ Premium V10

Blood Glucose Monitoring System

Owner's Manual

Dear ForaCare Premium V10 System Owner:

Thank you for purchasing the **ForaCare Premium V10** Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read and follow the instructions in this Owner's Manual and any package inserts for test strips and control solution.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the **ForaCare Premium V10** Blood Glucose Monitoring System to easily monitor your blood glucose levels.

If you have questions, just call ForaCare Customer Care Service Center toll-free at 1-866-469-2632. We are here for you 24 hours a day, 7 days a week. You can also visit www.foracare.com.

IMPORTANT SAFETY PRECAUTIONS

READ BEFORE USE

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

1. Use this device **ONLY** for the intended use described.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
5. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
6. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
7. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

BEFORE YOU BEGIN	06
Important Information	06
Intended Use	07
Test Principle	07
Meter Overview	08
Display Screen	09
Test Strip	10
Setting the Meter	11
THE TWO MEASURING MODES	14
Control Solution Testing	15
Performing a Control Solution Test	15
TESTING WITH BLOOD SAMPLE	18
Preparing the Puncture Site	18
Alternative Site Testing	20
Performing a Blood Glucose Test	21
METER MEMORY	24
Reviewing Test Results	24
Reviewing Blood Glucose Day Average Results	25
Downloading Results onto a Computer	26
MAINTENANCE	27
Battery	27
Caring for Your Meter	28
Caring for Your Test Strips	31
Important Control Solution Information	32
SYSTEM TROUBLESHOOTING	33
Result Readings	33
Error Messages	34
Troubleshooting	36
DETAILED INFORMATION	37
SPECIFICATIONS	38

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than normal, and you do not have any symptoms of illness, repeat the test. If you have symptoms or continue to get results which are higher or lower than normal, consult with your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings lower than actual values may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult a healthcare professional before use.
- Critically ill patients should not be tested with a blood glucose monitoring system.

Intended Use

The **ForaCare Premium V10** Blood Glucose Monitoring System is intended for use in the quantitative measurement of glucose in fresh capillary whole blood from the finger and the following alternative sites: palm, forearm, upper-arm, calf and thigh. It is intended for use by people with diabetes mellitus at home as an aid in monitoring the effectiveness of diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and is not intended for use on neonates. Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing.


The alternative site testing in the **ForaCare Premium V10** Blood Glucose Monitoring System can be used only during steady-state blood glucose conditions.

This meter has some speaking functions but is not intended for use by the visually impaired.

Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Meter Overview

1 Test Strip Slot

Insert test strip here to turn the meter on for testing.

2 Test Strip Ejector

Eject the used strip by sliding up this button.

3 Display Screen

4 M Button

Display the meter memory; silence a reminder alarm, change a setting in SET mode.


5 Data Port

Download test results with a cable connection.

6 SET Button


Enter SET mode and confirm the meter settings.

Display Screen

- | | |
|----------------------------|---------------------------------|
| 1 Volume Symbol | 8 Memory Symbol |
| 2 Insert Test strip | 9 Ketone Warning |
| 3 Test Strip Symbol | 10 Day Average |
| 4 Test Result | 11 Control Solution Mode |
| 5 Blood Drop Symbol | 12 Measurement Unit |
| 6 Date | 13 Low Battery Symbol |
| 7 Time | |

Test Strip

ATTENTION:

The front side of test strip should face up when inserting test strip.


Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE:

The **ForaCare Premium V10** monitor should only be used with **ForaCare Premium V10** Test Strips. Using other test strips with this meter can produce inaccurate results.


Setting the Meter

Before using your meter for the first time or if you change the meter battery, you should check and update the settings. Make sure you complete the steps below and have your desired settings saved.

Entering the Setting Mode

Start with the meter off (no test strip inserted). **Press SET.**


1. Setting the date

With the year flashing, press **M** until the correct year appears. Press SET.

With the month flashing, press **M** until the correct month appears. Press SET.

With the day flashing, press **M** until the correct day appears. Press SET.

2. Setting the time



With the hour flashing, press **M** until the correct hour appears. Press SET.

With the minute flashing, press **M** until the correct minute appears. Press SET.

3. Deleting the memory

With “dEL” and a flashing “M” on the display, press **M**. If recorded results does not need to be deleted from the memory, press SET to skip.

With “dEL” and a flashing “M” on the display, press **M** to delete all readings in memory.



Screen will display “—” and “☺”, press SET to set the speaking volume, press **M** to turn-off the meter.

4. Setting the speaking volume

There are eight (8) speaking volume levels to choose from. Press **M** until the desired speaking volume appears. To confirm your selection, press SET.

Volume 0 indicates that the speaking function is turned off, and “🔊” will not display during testing. Volume 1 to 7 indicates speaking volume from low to high, and “🔊” will be displayed throughout the testing.

5. Choosing a language

Press **M** to select L1/L2. The default language for the meter is L1, which is English. L2 is the alternate language as indicated in your packaging. To confirm your selection, press SET.

Congratulations! You have completed all settings!

NOTE:


- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idles for 3 minutes during the setting mode, it will switch off automatically.

THE TWO MEASURING MODES


The meter provides you with two modes for measuring, General and QC.

MODES	USE WHEN
General	any time of day without regard to time since last meal
QC	testing with the control solution

You can switch between the two modes mode by:

1. With the meter turned off, insert a test strip to turn on the meter. The screen will display a flashing blood drop.

2. Press **M** to switch between General and QC mode.

Control Solution Testing


Our control solutions contain a known amount of glucose that reacts with test strips and is used to ensure that the meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have all items prior to testing your blood glucose.

Test with a control solution when:

- you first receive the meter,
- you first receive the meter.
- you do your weekly routine quality testing of the meter and test strips.
- you begin using a new vial of test strips.
- you suspect the meter and or test strips are not working correctly.
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate.
- you are practicing the testing procedure.
- you have dropped or think you may have damaged the meter.

Performing a Control Solution Test


1. Insert the test strip to turn on the meter

Wait for the meter to display “” and “”.

Voice Guide

Thank you for using this product. Please relax during measurement. Please apply blood into the test strip.


2. Press **M** until “CTL” is displayed.

With “CTL” displayed, the meter will not store your test result in memory. If you press **M** again, the “CTL” will disappear and this test is no longer a control solution test.

WARNING:

When doing the control solution test “CTL” must be displayed so that the test results will not be stored in the memory. Failure to have the meter in the correct mode will result in having the control solution test results stored with your actual blood glucose test results. This will cause an inaccurate analysis of your stored results.


3. Apply control solution

Shake the control solution bottle thoroughly before use. Squeeze out a drop onto the tip of the bottle cap and wipe it off using a clean tissue or cloth. Then squeeze out another drop onto the tip of the bottle cap.

Hold the meter so that the strip is pointing downward and move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down.

To avoid contaminating the control solution, do not directly apply control solution onto a strip or otherwise touch the nozzle of the bottle.

4. Read and compare the result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial. It should fall within the range. If not, please read the instructions again and repeat the control solution test.

Voice Guide

The blood glucose is (number) milligrams per deciliter.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood using this meter. Go to www.foracare.com, or call your customer service for help.

NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **Meter Maintenance** section for important information about your control solutions.

TESTING WITH BLOOD SAMPLE

Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets or the lancing device.
- Wash and dry your hands thoroughly before and after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by that person.

Sharing the lancing device and lancets may increase the risk of spreading infectious diseases. Lancing device must not be used on more than one person.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from an alternate site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:


Wash and dry your hands before starting.

- Select the puncture site either at fingertips or another body part (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap while setting up the lancing device.

• **Fingertip testing**

Press the tip of the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger, a click indicates that the puncture is complete.


• **Blood from sites other than the fingertip**

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm, hand, thigh, or calf, avoid lancing the areas with obvious veins to avert excessive bleeding.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing

Important:

There are limitations with Alternative Site Testing (AST).

Please consult your health care professional before you perform AST.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, **take a blood sample from your finger only.**


We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not match the way you feel.
- Your routine glucose results often fluctuate.

Performing a Blood Glucose Test

1. Insert the test strip to turn on the meter

Wait for the meter to display “” and flashing “”.


Voice Guide

Thank you for using this product. Please relax during measurement. Please apply blood into the strip.

2. Select the appropriate measuring mode (General or QC) by pressing .

For selecting the measurement mode, please refer to the “TWO MEASURING MODES”.


3. Obtaining a blood sample

Use the pre-set lancing device to puncture the chosen site. After lancing wipe away the first drop of blood with a clean cotton swab. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

The volume of blood sample must be at least **0.5** microliter (μL) of volume.

( actual size)

4. Apply the sample

Hold the meter with the strip pointing downward and the absorbent hole of the test strip touching the blood drop. Blood will be drawn in to the strip. Once the confirmation window is completely filled, meter begins counting down. Do not remove the strip from the blood drop until the meter beeps.

NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or our customer care service for assistance.


5. Read your result

The result of your blood glucose test will appear after the meter counts down to 0. This blood glucose result will automatically be stored in the memory.

Voice Guide

The blood glucose is (number) milligrams per deciliter.

6. Eject the used test strip and remove the lancet

To eject the test strip, point the strip into a disposal container for contaminated items. Slide the Test Strip Ejector forward. The meter will automatically turn itself off after the test strip is ejected.

When removing the lancet, follow the instructions from the lancing device insert.


WARNING:

The used lancet and test strip may be biohazardous. Please follow the proper procedure for discarding biohazardous waste in your area.

METER MEMORY

The meter stores up to 450 most recent blood glucose test results with date and time stamp in its memory. To enter the meter memory function, **start with the meter turned off.**

Reviewing Test Results



1. Press and release **M** “**M**” will appear on the display. Press **M** again, and the first reading you see is the last blood glucose result along with date, time and the measurement mode.

Voice Guide

The blood glucose is (number) milligrams per deciliter.


2. Press **M** **again** to recall the next test results stored in the meter each time you press. **After the last test result, press** **M** **again and the meter will turn off.**

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
Reviewing Blood Glucose Day Average Results

1. Press and release **M.** When “**M**” appears on the display, keep pressing **M** for 3 seconds until the flashing “**DAY AVG**” appears. Release **M** and then your 7-day average result measured in general mode will appear on the display.

2. Press **M to review** 14-, 21-, 28-, 60- and 90- day average results stored in memory.

3. Exit the meter memory. Press **M** after the 90-day average is displayed and the meter will turn off.

NOTE:

- Any time you wish to exit the memory, hold press **M** down for 5 seconds or leave the meter idle for 3 minutes. The meter will turn off automatically.
- Control solution results are **NOT** included in the day average.
- If there are no records in memory, “—” displays when you recall the test results or review the average result.


Downloading Results onto a Computer

Data Transmission Via Cable

You can use the meter with the ForaCare 4-pin USB cable and the Health Care System Software to view your test results on your personal computer. To learn more about the Health Care System Software or to obtain an USB cable, please visit www.foracare.com, or call our customer care service number for further assistance.

1. Obtaining the required cable and installing the software

To download the Health Care System Software, please visit the ForaCare Inc. website at www.foracare.com.

2. Connecting to a personal computer

Connect the cable to a USB port on your computer. With the meter turned off, connect the other end of the USB cable to the meter data port. “USB” will appear on the meter display, indicating that the meter is in communication mode.

3. Data transmission

To transmit data, follow the instructions provided with the software. Test results will be transmitted with date and time stamp. After the results have been transmitted, remove the cable, meter will automatically turn off.

WARNING:

While the meter is connected your computer, it will not be able to perform a blood glucose test.


MAINTENANCE

Battery

Your meter comes with two 1.5V AAA size alkaline batteries.

Low Battery Signal

The meter will display the message below to alert you when the meter power is getting low.


The “” appears with E-b and low:

The power is not enough to perform testing. Please change the batteries immediately.

Replacing the Battery

To replace the batteries, make sure that the meter is turned off.

1. Remove the battery cover.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Close the battery cover, paying close attention to the polarity of the batteries. (If the batteries are inserted correctly, you will hear a “beep”).

NOTE:

- Replacing the batteries does not affect the test results stored in the memory.

- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Caring for Your Meter

Wash and dry your hands thoroughly before use or handling of the meter and strips.

What is Cleaning and Disinfection?

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter?

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter?

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. **Do NOT use organic solvents to clean the meter.**

We recommend using the product listed below to clean and disinfect your meter. The recommended product contains isopropyl alcohol as its active

ingredient. It has been shown to be safe when used on the ForaCare Premium V10 blood glucose meter.

- Micro-Kill+™ (Micro-Kill Plus™) by Medline
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all exterior surface of the meter including the display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot, be careful not to allow excess liquid to get inside. Keep the meter moist with disinfection solution contained in the wipe for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe. Use two or more wipes if necessary.

3. Remove the wipe. Allow the meter surface to dry completely.

4. Discard the used wipes.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at

1-866-469-2632 for assistance.

This device has been validated to withstand 5,000 disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 1,000 tests per year on each meter and around 3 disinfection cycles per meter per day for 5 years, the expected life of device.

Stop using the meter if you see any signs of deterioration such as the LCD display is cracked or cloudy, the buttons no longer function, or the outer casing is cracked. Please contact the customer service at 1-866-469-2632 for a replacement meter if any of the signs of deterioration are noticed.

The meter should be replaced after the validated number of disinfection cycles or the warranty period, whichever comes first.

NOTE:

- Do **NOT** clean or disinfect the meter while performing a test.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to handling by an assistant.
- Do **NOT** allow cleaning or disinfecting solution to get in the test slot, battery compartment, or strip-ejection button.
- If you do get moisture in the test strip slot, wipe it away with the corner of a tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case and damage the electronic components or circuitry.

Meter Storage

- Storage conditions: -4 °F to 140 °F (-20 °C to 60 °C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage conditions: 39.2 °F to 104 °F (4 °C to 40 °C), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months of first opening.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution bottle and discard the remaining solution after 3 months of first opening.
- It is recommended that the control solution test be done at room temperature (68°F to 77°F / 20°C to 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the bottle before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6°F and 86°F (2°C and 30°C). Do **NOT** freeze.

SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair the meter yourself and never try to disassemble the meter under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
 The image shows a digital display with 'Lo' at the top, a small speaker icon in the top right, '8-10' in the middle, and '10:00 AM' at the bottom.	Appears when your blood glucose test result is below the lower measurement limit, which is less than 20 mg/dL.
 The image shows a digital display with '388' at the top, a small speaker icon in the top right, 'mg/dL' below the number, '8-10' in the middle, and '10:00 AM' at the bottom.	Appears when your result is equal to or higher than 240 mg/dL. This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
 The image shows a digital display with 'Hi' at the top, a small speaker icon in the top right, '8-10' in the middle, and '10:00 AM' at the bottom.	Appears when your blood glucose test result is higher than the upper limit of measurement, which is 600 mg/dL.

Error Messages

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	WHAT TO DO
	Only Error Beep	Problem with the meter.	Review the instructions and repeat the test. If the meter still does not work, contact the local distributor for assistance.
	Only Error Beep		
	Only Error Beep		
	Only Error Beep		
	Only Error Beep	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.
	Only Error Beep	Appears when the batteries cannot provide enough power for a test.	Replace the batteries immediately.

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	WHAT TO DO
	Only Error Beep	Appears when test strip is removed while counting down, or insufficient blood volume.	Repeat the test with a new test strip.
	Only Error Beep	Appears when a used test strip is inserted.	Repeat with a new test strip.
	Only Error Beep	Appears when ambient temperature is below system operation range.	System operation range is 50 °F to 104 °F (10 °C to 40 °C). Repeat the test after the meter and test strip are in the above temperature range.
	Only Error Beep	Appears when ambient temperature is above system operation range.	

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer service.

2. If the test does not start counting down after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (20°C to 25°C/68°F to 77°F) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer service.

DETAILED INFORMATION

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2010). Clinical Practice Recommendations. Diabetes Care, 33 (Supplement 1): S1–S100.

Please consult your doctor to determine a target range that works best for you.

SPECIFICATIONS

Model No.: TD-4124

Dimension & Weight: 88.67 x 55.61 x 23.07 mm, 48.8 g

Power Source: Two 1.5V AAA alkaline batteries

Display: LCD

Memory: 450 measurement results with respective date and time

External output: USB

Auto electrode insertion detection

Auto sample loading detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature Warning

Operating Condition: 50 °F to 104 °F (10 °C to 40 °C), below 85% R.H.
(non-condensing)

Storage/Transportation Conditions: -4 °F to 140 °F (-20 °C to 60 °C),
below 95% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

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1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)

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Fora Care Inc.